



PM-7

Features

- 7 inch high-resolution color screen.
- Exquisite design, suitable for moving operation.
- Various interfaces: standard screen, trend screen, oxyCRG screen, big front screen.
- Real-time S-T segment analysis, pace-maker detection.
- Up to 400 groups NIBP list, 60 alarm event records recall, 7-day trend chart in storage.
- Indication of incorrect manipulation and malfunction analysis.
- Capture 2-hour dynamic trend.
- Built-in rechargeable lithium battery, 2-hour working capability.

Standard Configuration

SpO2, PR, NIBP, ECG, RESP, TEMP

Optional Configuration

ETCO2



PM-8B

Features

- 8 inch high-resolution color screen.
- Exquisite design, suitable for moving operation.
- Various interfaces: standard screen, trend screen, oxyCRG screen, NIBP list screen, big front screen.
- Real-time S-T segment analysis, pace-maker detection.
- Up to 400 groups NIBP list, 60 alarm event records recall in storage.
- Built-in rechargeable lithium battery, 2-hour working capability.
- Capture 2-hour dynamic trend.

Standard Configuration

SpO2, PR, NIBP, ECG, RESP, TEMP

Optional Configuration

ETCO2

